

Exchange Guide Wire acc. to Schueller

Modified exchange guide-wire

This guide-wire offers superior stability in comparison to conventional Seldinger guide-wires. Pigtail drainage catheters can be easily straightened and introduced. In particular by catheter exchange, this guide-wire proves to be an essential instrument for the urologist. The flexible end facilitates the introduction of various instruments.

Available in straight or J-curved form.

Exchange guide wire straight, flexible end, rigid shaft (40 cm), flexible tip (10 cm)



Diameter	Length	Order-Nr.
inch	cm	
.035	100	1135-1100

Diameter	Length	Order-Nr.
inch	cm	
.035	90	1135-1290
.035	100	1135-1300