

BigLumen Aspiration Catheter

for the percutaneous Aspiration-Thromboembolectomy to ease the removal of thromboembolic material

BigLumen - Aspiration catheter, straight

perfect supplement

Percutaneous Introducer Sheaths with haemostatic valve



Order-Nr.	Length		Inner Diameter
	cm	F	
1720-0100	100	5	1.25

1720-0120	120	5	1.25
1720-1070	70	6	1.43
1720-1100	100	6	1.43
1720-2090	90	7	1.64
1720-2100	100	7	1.64
1720-3080	80	8	1.88
1720-3100	100	8	1.88

Box: 5 PCS