


TIPS Puncture Set „Special I“

- Ideal torque stability
- 30° or 60° curved needle configuration
- Ultrasound visibility due to sandblasting
- Maximized sharpness of needle due to Colapinto-cut
- Perfect adaption of guiding catheter, puncture needle and accessory
- Optimal configuration and stability of guiding catheter



Order Code	Set Description Version „Special 1“		Order Code single items	
1490-6000	Needle 1.8 mm with Obturator Type I - 30 ° curved	Box: 1 unit	1490-6100	Box: 1 unit
	Special guide wire .035 inch - 180 cm - straight		1490-1500	Box: 5 units
	8F Guiding catheter , 50 cm, adapted to needle		1490-1200	Box: 5 units
1490-7000	Needle 1.8 mm with Obturator Type II - 60 ° curved	Box: 1 unit	1490-7100	Box: 1 unit
	Special guide wire .035 inch - 180 cm - straight		1490-1500	Box: 5 units
	8F Guiding catheter , 50 cm, adapted to needle		1490-1200	Box: 5 units