


AltaFlow Measuring Catheter

for intravascular measurements
 can be used as a flush-catheter at the same time
 up to a burst pressure of 1,050 PSI
 high radiopacity
 tungsten markers and tip for use in MR procedures



	Length	Side Holes	Measuring Length	Adapted to GW	Flow Rate	Pressure	Order-Nr.
F	cm		cm	inch	ml/s	PSI	
5	90	8	14	.035	13	1050	1080-1000
5	90	8	22	.035	13	1050	1080-2000

	Length	Side Holes	Measuring Length	Adapted to GW	Flow Rate	Pressure	Order-Nr.
F	cm		cm	inch	ml/s	PSI	
4	110	8	22	.035	12	1050	1080-0200
5	90	8	14	.035	13	1050	1080-3000
5	90	8	22	.035	13	1050	1080-4000
5	110	8	22	.035	12	1050	1080-4100