


## Biopsy Needle

**Biopsy needle compatible with the Vitesse biopsy system.**

Cannula with sand-blasted tip for better ultrasound visibility.  
Centimetre markers facilitate CT guided procedures.  
Locking device for fixing the stylet

Biopsy needle





Biopsy-Needle				Coaxial Needle / Trocar, 2-part			
Length	Diameter		Order code	Length	Diameter		Order code
mm	mm	Gauge		mm	mm	Gauge	
100	1,2	18	1399-1210	-	-	-	-
100	1,6	16	1399-1610	-	-	-	-
100	2,1	14	1399-2110	70,5	2,4	13	1397-2110
150	1,2	18	1399-1215	120,5	1,5	17	1397-1215
150	1,6	16	1399-1615	-	-	-	-
200	1,2	18	1399-1220	170,5	1,5	17	1397-1220
200	1,6	16	1399-1620	-	-	-	-
280	1,2	18	1399-1228	-	-	-	-
<b>Box: 10 units</b>				<b>Box: 10 units</b>			