




Ureteral Stent Set (OptiSoft)

- maximum patient comfort due to extra-soft material
- for short and long term indication
- smooth surface with low incrustation tendency
- high x-ray visibility
- excellent drainage performance (large lumen, extra-large side holes in the pigtail and/or in the shaft)

Available in the versions:

- ureteral stent set, both end open, with guide wire
- ureteral stent set, both end open, without guide wire
- ureteral stent set, distal end closed-Tiemann, with guide wire





ureteral stent set, both end open, with guide wire

Catheter Ø	Length A	Length B	GW-Ø	length guid wire	Order-Nr.
CH	cm	cm	inch	cm	
4,5	26	22	.032	150	3004-0300
4,5	28	24	.032	150	3004-0400
4,5	30	26	.032	150	3004-0500
4,5	32	28	.032	150	3004-0600*
6	26	22	.035	150	3004-1300
6	28	24	.035	150	3004-1400
6	30	26	.035	150	3004-1500
6	32	28	.035	150	3004-1600
6	34	30	.035	150	3004-1700*
6	36	32	.035	150	3004-1800*
7	26	22	.035	150	3004-2300
7	28	24	.035	150	3004-2400
7	30	26	.035	150	3004-2500
7	32	28	.035	150	3004-2600
7	34	30	.035	150	3004-2700*
7	36	32	.035	150	3004-2800*
8	26	22	.038	150	3004-3300*
8	28	24	.038	150	3004-3400*
8	30	26	.038	150	3004-3500
8	34	30	.038	150	3004-3700*
8	36	32	.038	150	3004-3800*

Box: 10 PCS

ureteral stent set, both end open, without guide wire

Catheter Ø	Length A	Length B	Adapted to guide wire	Order-Nr.
CH	cm	cm	inch	
4,5	26	22	.032	3004-0330*
4,5	28	24	.032	3004-0430
4,5	30	26	.032	3004-0530
4,5	32	28	.032	3004-0630
6	26	22	.035	3004-1330
6	28	24	.035	3004-1430
6	30	26	.035	3004-1530*
6	32	28	.035	3004-1630
6	34	30	.035	3004-1730*
6	36	32	.035	3004-1830*
7	26	22	.035	3004-2330*
7	28	24	.035	3004-2430
7	30	26	.035	3004-2530
7	34	30	.035	3004-2730*
7	36	32	.035	3004-2830*
8	26	22	.038	3004-3330*
8	28	24	.038	3004-3430*
8	30	26	.038	3004-3530*

Box: 10 PCS

ureteral stent set, distal end closed-Tiemann, with guide wire

Catheter Ø	Length A	Length B	GW-Ø	length guid wire	Order-Nr.
CH	cm	cm	inch	cm	

Catheter Ø	Length A	Length B	GW-Ø	length guid wire	Order-Nr.
CH	cm	cm	inch	cm	
4,5	26	22	.032	150	3003-0300
4,5	28	24	.032	150	3003-0400
4,5	30	26	.032	150	3003-0500*
4,5	32	28	.032	150	3003-0600*
6	26	22	.035	150	3003-1300*
6	28	24	.035	150	3003-1400*
6	30	26	.035	150	3003-1500*
6	32	28	.035	150	3003-1600*
6	34	30	.035	150	3003-1700*
6	36	32	.035	150	3003-1800*
7	26	22	.035	150	3003-2300*
7	28	24	.035	150	3003-2400*
7	30	26	.035	150	3003-2500*
7	32	28	.035	150	3003-2600
7	34	30	.035	150	3003-2700*
7	36	32	.035	150	3003-2800*
8	26	22	.038	150	3003-3300*
8	28	24	.038	150	3003-3400*
8	30	26	.038	150	3003-3500*
8	32	28	.038	150	3003-3600*
8	34	30	.038	150	3003-3700*
8	36	32	.038	150	3003-3800*

Box: 10 PCS