


Pediatric Ureteral Stent Set (OptiSoft)

Complete set consisting of: ureteral stent, guide-wire and pusher

Extra short ureteral stent for pediatrics. Ureteral stent diameter in 3 CH and 4 CH.
 Guide wire and pusher are adapted to the ureteral stent diameter.



ureteral stent set, both end open, with guide wire

Catheter Ø	Length A	Length B	GW-Ø	length guid wire	Order-Nr.
3CH	30cm	30cm	1.2mm		

Catheter Ø	Length A	Length B	GW-Ø	length guid wire	Order-Nr.
CH	cm	cm	inch	cm	
3	10	8	.021	100	3001-1010
3	12	10	.021	100	3001-1020
3	14	12	.021	100	3001-1030
3	16	14	.021	100	3001-1040
3	18	16	.021	100	3001-1050*
3	20	18	.021	100	3001-1060*
3	22	20	.021	100	3001-1070*
4	10	8	.025	100	3001-1100*
4	12	10	.025	100	3001-1200
4	14	12	.025	100	3001-1300*
4	16	14	.025	100	3001-1400

Box: 5 PCS

ureteral stent set, distal end closed-Tiemann, with guide wire

Catheter Ø	Length A	Length B	GW-Ø	length guid wire	Order-Nr.
CH	cm	cm	inch	cm	
3	10	8	.021	100	3001-0010*
3	12	10	.021	100	3001-0020*
3	14	12	.021	100	3001-0030*
3	16	14	.021	100	3001-0040*
4	10	8	.025	100	3001-0100*
4	12	10	.025	100	3001-0200*
4	14	12	.025	100	3001-0300*
4	16	14	.025	100	3001-0400*

Box: 5 PCS