




## Steerable ureteral stent set, optimatic, multilength (OptiSoft)

Steerable ureteral stent with double pigtail for variable lengths and simple application system - facilitates the manipulation (turning, pushing, pulling) within the urinary tract

*assembled system*





ureteral stent set, both end open, with guide wire

Catheter Ø	Length A	Length B	GW-Ø	length guid wire	Order-Nr.
CH	cm	cm	inch	cm	
6	24-32	20-28	.035	150	3014-2060
7	24-32	20-28	.035	150	3014-3060

Catheter Ø	Length A	Length B	GW-Ø	length guid wire	Order-Nr.
CH	cm	cm	inch	cm	
8	24-32	20-28	.038	150	3014-4060*

Box: 10 PCS

ureteral stent set, both end open, without guide wire

Catheter Ø	Length A	Length B	Adapted to guide wire	Order-Nr.
CH	cm	cm	inch	
6	24-32	20-28	.035	3014-2065*
7	24-32	20-28	.035	3014-3065*
8	24-32	20-28	.038	3014-4065*

Box: 10 PCS