



Universal Nitinol Stone Basket, helical

Universal stone basket for removal of urological stones.

The nitinol stone basket (helical with 4 wires) consists of three primary components: basket, sheath and handgrip.

Advantages:

- Removable handgrip / sheath
- Single-handed operation
- Excellent basket expansion
- Torqueable 1:1
- Excellent pushability



Sheath Ø	Sheath Length	Ø Basket	Basket length	Order-Nr.
CH	cm	mm	mm	
2,7	90	12	25	3500-0210
4	90	15	35	3500-1200

Box: 1 PCS