


optimed

**Opti-J**

Ureteral stent system

**2 in 1**


One ureteral  
stent with two  
possibilities of use



STEERABLE URETERAL STENT & URETERAL STENT WITH PUSHER

**Opti-J**





| CH  | Länge<br>longueur<br>length<br>- cm - |    | FD<br>FG<br>GW<br>Ø / inch | <br>Ureterschiene<br>mit Führungsdraht<br><i>stent urétéral<br/>avec fil-guide</i><br>ureteral stent<br>with guide wire | <br>Ureterschiene<br>ohne Führungsdraht<br><i>stent urétéral<br/>sans fil-guide</i><br>ureteral stent<br>without guide wire | <br>Ureterschiene<br>mit Führungsdraht<br><i>stent urétéral<br/>avec fil-guide</i><br>ureteral stent<br>with guide wire |
|---|---------------------------------------|----|----------------------------|---|---|---|
|   | A                                     | B  |                            | Set   | Mini Set  | Set   |
| 4,5   | 26                                    | 22 | .032                       | 8200-0600   | 8200-0630   | 8100-0600   |
|   | 28                                    | 24 |                            | 8200-0700   | 8200-0730   | 8100-0700   |
|   | 30                                    | 26 |                            | 8200-0800   | 8200-0830   | 8100-0800   |
|   | 32                                    | 28 |                            | 8200-0900   | 8200-0930   | 8100-0900   |
|   | 34                                    | 30 |                            | 8200-0950   | 8200-0940   | 8100-0950   |
| Länge Führungsdraht / longueur du fil-guide / length guide wire |                                       |    |                            | 150 cm  | --  | 125 cm  |
| 6   | 26                                    | 22 | .035                       | 8200-1500   | 8200-1530   | 8100-1500   |
|   | 28                                    | 24 |                            | 8200-1600   | 8200-1630   | 8100-1600   |
|   | 30                                    | 26 |                            | 8200-1700   | 8200-1730   | 8100-1700   |
|   | 32                                    | 28 |                            | 8200-1800   | 8200-1830   | 8100-1800   |
|   | 34                                    | 30 |                            | 8200-1850   | 8200-1840   | 8100-1850   |
|   | 36                                    | 32 |                            | 8200-1890   | 8200-1920   | 8100-1860   |
| 7   | 26                                    | 22 | .035                       | 8200-1900   | 8200-1930   | 8100-1900   |
|   | 28                                    | 24 |                            | 8200-2000   | 8200-2030   | 8100-2000   |
|   | 30                                    | 26 |                            | 8200-2100   | 8200-2130   | 8100-2100   |
|   | 32                                    | 28 |                            | 8200-2200   | 8200-2230   | 8100-2200   |
|   | 34                                    | 30 |                            | 8200-2250   | 8200-2240   | 8100-2250   |
|   | 36                                    | 32 |                            | 8200-2290   |   | 8100-2260   |
| 8   | 26                                    | 22 | .038                       | 8200-2300   | 8200-2330   | 8100-2300   |
|   | 28                                    | 24 |                            | 8200-2400   | 8200-2430   | 8100-2400   |
|   | 30                                    | 26 |                            | 8200-2500   | 8200-2530   | 8100-2500   |
|   | 32                                    | 28 |                            | 8200-2600   | 8200-2630   | 8100-2600   |
|   | 34                                    | 30 |                            | 8200-2650   | 8200-2640   | 8100-2650   |
|   | 36                                    | 32 |                            | 8200-2690   |   | 8100-2660   |
| Länge Führungsdraht / longueur du fil-guide / length guide wire |                                       |    |                            | 150 cm  | --  | 150 cm  |



## 2 in 1

Eine Ureterschiene mit zwei Anwendungsmöglichkeiten:

- Als steuerbare Ureterschiene
- Als Ureterschiene mit Schieber

## 2 en 1

Un stent urétéral avec deux possibilités d'utilisation:

- Comme stent urétéral connectable et orientable
- Comme stent urétéral avec poussoir

## 2 in 1

One ureteral stent with two possibilities of use:

- Steerable ureteral stent
- Ureteral stent with pusher



### CLOSED DISTAL END WITH TUOHY-BORST ADAPTER



Ureterschiene, distales Ende geschlossen mit Tuohy-Borst-Adapter zur Fixierung des PTFE-Führungsdrahtes mit beweglicher Seele

Stent urétéral, extrémité distale fermée avec adaptateur Tuohy-Borst pour fixation au guide à âme mobile revêtu de PTFE

Ureteral stent, distal end closed with Tuohy Borst Adapter for fixation of PTFE coated guide wire with movable core

Verpackungseinheit: 5 Stück  
conditionnement: 5 unités  
box: 5 units